

Anti-Phospho-PTPN1-Tyr-152 antibody (STJA0005554)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Antigen Affinity Purified |
| Dilution Range | WB 1:1000 |
| Formulation | PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|--|
| Gene ID | 5770 |
| Gene Symbol | PTPN1 |
| UniProt ID | PTN1_HUMAN |
| Specificity | This antibody was cross-adsorbed to a phospho-tyrosine peptide then affinity purified with phospho-PTP1B (Tyr-152) peptide. The purified antibody detects a 50 kDa* band corresponding to PTP1B in western blots of pervanadate treated human A431 and Jur |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|