

## Anti-BRAF antibody (N-Term) (STJA0005558)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">673</a>
Gene Symbol	<a href="#">BRAF</a>
UniProt ID	<a href="#">BRAF_HUMAN</a>
Immunogen Region	N-Term
Specificity	The antibody detects 105 and 95 kDa* bands corresponding to the molecular mass of B-Raf. The doublet observed in Jurkat cells may represent multiple isoforms of the B-Raf, since all bands are blocked specifically with B-Raf (N-terminus) peptide.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---