

Anti-Phospho-BRAF-Ser-446 antibody (STJA0005559)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	673
Gene Symbol	BRAF
UniProt ID	BRAF_HUMAN
Specificity	This antibody was cross-adsorbed to phospho-B-Raf (Ser-579) peptide before affinity purification using phospho-B-Raf (Ser-446) peptide. The purified antibody detects a strong band at 95 kDa* and a weaker band at 105 kDa in western blots of human Jurk

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---