

Anti-RCAN1 antibody (central Region) (STJA0005563)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 6:00
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1827
Gene Symbol	RCAN1
UniProt ID	RCAN1_HUMAN
Immunogen Region	central Region
Specificity	The antibody detects several isoforms of RCAN1 on SDS-PAGE immunoblots of mouse C2C12, rat PC12, and human Jurkat cells. The predominate isoform detected is the 38 kDa RCAN1L, but the antibody may also detect RCAN2 and RCAN3.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---