

Anti-SEMA6A antibody (C-Term) (STJA0005577)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	57556
Gene Symbol	SEMA6A
UniProt ID	SEM6A_HUMAN
Immunogen Region	C-Term
Specificity	The antibody detects a 120kDa* protein corresponding to the apparent molecular mass of Sema-6A on SDS-PAGE immunoblots of Sema-6A transfected COS7 cells, human colorectal cancer cells, and rat PC12 cells. In addition, shRNA knockdown of Sema-6A expre
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.