

Anti-SHC1 antibody (C-Term) (STJA0005579)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6464
Gene Symbol	SHC1
UniProt ID	SHC1_HUMAN
Immunogen Region	C-Term
Specificity	The antibody detects 46/52/66 kDa* bands corresponding to the molecular weight of ShcA variants on SDS-PAGE immunoblots of A431 cells. It also detects ShcA in human Jurkat, HeLa, and SKN-SH cells, as well as rat brain tissue and rabbit spleen fibrobl

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---