

Anti-Phospho-SSH1-Ser-978 antibody (STJA0005585)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54434
Gene Symbol	SSH1
UniProt ID	SSH1_HUMAN
Specificity	The antibody detects a 150 kDa* full length recombinant human SSH1L phosphoprotein and this reactivity is eliminated with lambda phosphatase treatment. The antibody also detects SSH1L on immunoblots of rat PC12 and mouse C2C12 cells that were stimulated.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---