

## Anti-SMPD4 antibody (C-Term) (STJA0005588)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:250 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">77626</a>
Gene Symbol	<a href="#">Smpd4</a>
UniProt ID	<a href="#">NSMA3_MOUSE</a>
Immunogen Region	C-Term
Specificity	The antibody detects an 85-90 kDa* band corresponding to nSMase3 on SDS-PAGE immunoblots of adult mouse heart and mouse C2C12 cells. The antibody works for western blot, ELISA, immunoprecipitation, and immunocytochemistry.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---