

Anti-Phospho-SOX2-Thr-118 antibody (STJA0005590)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ICC
	Host / Source	Rabbit
	Reactivity	Chicken/Human/Mouse/Rat/Xenopus
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Conjugation	Unconjugated
	Purification	Antigen Affinity Purified
	Dilution Range	WB 1:1000 ICC 1:50
	Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
	Isotype	IgG
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	20674
	Gene Symbol	Sox2
	UniProt ID	SOX2_MOUSE
	Specificity	The antibody was cross adsorbed to unphosphorylated Sox2 (Thr-118) peptide before affinity purification using phospho-Sox2 (Thr-118) peptide. This antibody detects a 34 kDa* protein on SDS-PAGE immunoblots of mouse F9 stem cells. This reacti
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.