

Anti-Methyl-SOX2-Lys-119 antibody (STJA0005591)

GENERAL INFORMATION		
Product Type	Primary antibodies	
Applications	WB/ICC	
Host / Source	Rabbit	
Reactivity	Chicken/Human/Mouse/Rat/Xenopus	
PRODUCT PROPERTIES		
Clonality	Polyclonal	
Conjugation	Unconjugated	
Purification	Antigen Affinity Purified	
Dilution Range	WB 1:1000 ICC 1:100	
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol	
Isotype	IgG	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.	
TARGET INFORMATION		
Gene ID	20674	
Gene Symbol	Sox2	
UniProt ID	SOX2_MOUSE	
Specificity	The antibody was cross adsorbed to unmethylated Sox2 (Lys-119) peptide before affinity purification using methylated-Sox2 (Lys-119) peptide. This antibody detects a 34 kDa* protein on SDS-PAGE immunoblots of mouse F9 stem cells, as well as h	
ADDITIONAL INFORMATION		
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.	