

Anti-TACSTD2 antibody (Extracellular region) [M005] (STJA0005615)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4070
Gene Symbol	TACSTD2
UniProt ID	TACD2_HUMAN
Immunogen Region	Extracellular region
Specificity	Clone M005 detects a 50-65 kDa* band corresponding to the molecular mass of Trop-2 on SDS-PAGE immunoblots of native human A431, MDA-MB-231, and MCF-7 cell lysates. The antibody also detects denatured Trop-2, but with less affinity. Lower dilutions o

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---