

## Anti-TRPM7 antibody (Extracellular region) (STJA0005617)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">54822</a>
Gene Symbol	<a href="#">TRPM7</a>
UniProt ID	<a href="#">TRPM7_HUMAN</a>
Immunogen Region	Extracellular region
Specificity	This antibody detects a 220 kDa* protein on SDS-PAGE immunoblots of human MCF7, MDA-MB-231, mouse brain, and rat PC12 cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---