

## Anti-TUBB3 antibody (central Region) (STJA0005629)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">10381</a>
Gene Symbol	<a href="#">TUBB3</a>
UniProt ID	<a href="#">TBB3_HUMAN</a>
Immunogen Region	central Region
Specificity	This antibody detects a 50 kDa* protein corresponding to the molecular mass of Beta-Tubulin on SDS-PAGE immunoblots of purified brain tubulin and mouse brain tissue.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---