

Anti-WASL antibody (STJA0005648)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8976
Gene Symbol	WASL
UniProt ID	WASL_HUMAN
Specificity	This antibody was cross-adsorbed to phospho-N-WASP (S484/S485) peptide (without carrier) then affinity purified using unphosphorylated N-WASP (S484/S485) peptide (without carrier). The antibody detects a 65 kDa* protein on SDS-PAGE immunoblots of rat

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---