

Anti-Phospho-WASF1-Tyr-125 antibody (STJA0005650)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8936
Gene Symbol	WASF1
UniProt ID	WASF1_HUMAN
Specificity	This antibody was cross-adsorbed to phospho-WAVE (Tyr-150) and unphosphorylated WAVE (Tyr-125) peptides before affinity purification using phospho-WAVE (Tyr-125) peptide. In western blots, the antibody detects an 80 kDa* band corresponding to phospho

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---