

Anti-NOS1 antibody (N-Term) (STJA0005670)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4842
Gene Symbol	NOS1
UniProt ID	NOS1_HUMAN
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of the human variant of nitric oxide synthase (NOS). The neuronal form of NOS (nNOS) is one of three forms and is enzymatically active as a dimer.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---