

## Anti-ARTN antibody (N-Term) (STJA0005674)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN <sub>3</sub>
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">9048</a>
Gene Symbol	<a href="#">ARTN</a>
UniProt ID	<a href="#">ARTN_HUMAN</a>
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of human artemin protein

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



antibody WB analysis Artemin immunolabeling in human dorsal root ganglia lysate STJA0005674