

Anti-IKBKG antibody (N-Term) (STJA0005690)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8517
Gene Symbol	IKBKG
UniProt ID	NEMO_HUMAN
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of human IKK gamma protein. Recognizes both human isoforms.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western Blot IKK-gamma immunolabeling in IL-1 stimulated 293 cell lysate STJA0005690