

Anti-TAFAZZIN antibody (STJA0005719)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6901
Gene Symbol	TAFAZZIN
UniProt ID	TAZ_HUMAN
Specificity	Reacts with the PDZ-binding domain of the 49 kDa human, mouse and rat TAZ protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot analysis MDCK cells demonstrating specific TAZ immunolabeling STJA0005719



Western Blot CTGF and YAP1-TAZ protein levels in 3D cancer spheroids STJA0005719