

Anti-RETNLA antibody (N-Term) (STJA0005735)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	57262
Gene Symbol	Retnla
UniProt ID	RETNA_MOUSE
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of mouse RELM alpha protein

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining RELM alpha green in mouse BV2 microglia cells with nuclear stain blue STJA0005735



Western Blot specific immunolabeling RELM alpha in transfected 293T cells STJA0005735



Western Blot specific immunolabeling RELM alpha in mouse BV2 cell lysate STJA0005735