

Anti-NLRC4 antibody (N-Term) (STJA0005739)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	58484
Gene Symbol	NLRC4
UniProt ID	NLRC4_HUMAN
Immunogen Region	N-Term
Specificity	Reacts with the N terminal region of the A and B isoforms of the human CLAN protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



antibody Western Blot CLAN-A and CLAN-B detection in lung lysate STJA0005739



Western Blot cN-II and Ipaf CARD12 protein expression in human cell lines STJA0005739