

Anti-TRPM8 antibody (C-Term) (STJA0005746)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	79054
Gene Symbol	TRPM8
UniProt ID	TRPM8_HUMAN
Immunogen Region	C-Term
Specificity	Reacts with the C-terminal region of human TRPMs protein

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining TRPM8 in formalin-fixed paraffin-embedded prostate cancer cells STJA0005746



Western Blot specific TRPM8 immunolabeling in transfected COS cells STJA0005746