

Anti-COL25A1 antibody (NC3 region) (STJA0005756)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	84570
Gene Symbol	COL25A1
UniProt ID	COPA1_HUMAN
Immunogen Region	NC3 region
Specificity	Reacts with the NC3 region of human and mouse CLAC-P.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining Collagen XXV in neocortex
Alzheimers diseased brain tissue STJA0005756



Western Blot specific immunolabeling Collagen
XXV alpha 1 in multiple lanes STJA0005756