

Anti-Shprh antibody (Internal Region) (STJA0005758)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268281
Gene Symbol	Shprh
UniProt ID	SHPRH_MOUSE
Immunogen Region	Internal Region
Specificity	Reacts with an internal region of the mouse SHPRH protein. 95% sequence identity to rat. 90% sequence identity to human

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western Blot analysis mouse brain lysate demonstrating SHPRH immunolabeling STJA0005758