

Anti-XRCC5 antibody (Internal Region) (STJA0005773)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Hamster/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN3
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7520
Gene Symbol	XRCC5
UniProt ID	XRCC5_HUMAN
Immunogen Region	Internal Region
Specificity	Reacts with the 80 kDa Ku80 protein. Sequence is 100% conserved between human, mouse, rat and hamster.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



immunoprecipitation IP Ku80 in CHO and A549 nuclear and cytoplasmic extracts STJA0005773



Western blot Ku80 protein 80 kDa in T47D and HeLa cells STJA0005773