

Anti-SLC2A10 antibody (Internal Region) (STJA0005777)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:250
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81031
Gene Symbol	SLC2A10
UniProt ID	GTR10_HUMAN
Immunogen Region	Internal Region
Specificity	Residues with internal residues of human, mouse and rat GLUT10 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining GLUT10 expression in mouse inner ear hair cells and spiral ganglion STJA0005777



Western Blot GLUT10 immunolabeling in liver and cochlea protein samples STJA0005777



Western Blot GLUT10 immunolabeling in transfected NCI-H1299 lung tumor cells STJA0005777