

Anti-UBIAD1 antibody (N-Term) (STJA0005778)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29914
Gene Symbol	UBIAD1
UniProt ID	UBIA1_HUMAN
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of the 37 kDa human TERE1 protein. This sequence is 85% identical in mouse and rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western Blot TERE1 immunolabeling in human bladder mucosa and urothelial carcinoma
 STJA0005778