

Anti-BCAS1 antibody (N-Term) (STJA0005783)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN3
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8537
Gene Symbol	BCAS1
UniProt ID	BCAS1_HUMAN
Immunogen Region	N-Term
Specificity	Reacts with the C-terminal region of human NABC1. NABC1 is a 60 kDa protein; however it forms dimers, which are easily detected in MCF-7 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



antibody Western Blot specific immunolabeling
NABC1 in MCF-7 cell lysate STJA0005783