

Anti-GRM2/3 antibody (C-Term) (STJA0005788)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:250
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2912
Gene Symbol	GRM2
UniProt ID	GRM2_HUMAN
Immunogen Region	C-Term
Specificity	Reacts with the C-terminal region of human, mouse and rat mGluR2/3 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot mGluR2-3 monomeric and dimeric forms in rat brain lysate STJA0005788