

## Anti-NR4A2 antibody (Internal Region) (STJA0005792)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN <sub>3</sub>
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">4929</a>
Gene Symbol	<a href="#">NR4A2</a>
UniProt ID	<a href="#">NR4A2_HUMAN</a>
Immunogen Region	Internal Region
Specificity	Reacts with an internal sequence of the human Nurr1 protein. Also reacts with rat and mouse Nurr1.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining Nurr1 in adult rat brain hippocampus tissue section STJA0005792



Western blot Nurr1 immunolabeling in human hippocampus tissue lysate STJA0005792