

## Anti-ADRA2C antibody (Internal Region) (STJA0005796)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC             |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Antigen Affinity Purified                                                                      |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:100-1:300                                                             |
| Formulation         | PBS + 0.02% NaN <sub>3</sub>                                                                   |
| Isotype             | IgG                                                                                            |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION


|                  |                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------|
| Gene ID          | <a href="#">152</a>                                                                              |
| Gene Symbol      | <a href="#">ADRA2C</a>                                                                           |
| UniProt ID       | <a href="#">ADA2C_HUMAN</a>                                                                      |
| Immunogen Region | Internal Region                                                                                  |
| Specificity      | Reacts with internal residues of the human Alpha 2c protein. Reacts with rat and mouse Alpha 2c. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining Alpha 2c Adrenergic Receptor in rat ventral spinal cord STJA0005796



Western Blot Alpha 2c Adrenergic Receptor in human brain lysate STJA0005796