

Anti-Tri-Methyl-H3-Lys37 antibody (STJA0005913)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	1 mg/mL
	Conjugation	Unconjugated
	Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
	Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
	Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
	Isotype	IgG
	Molecular Weight	Observed: 17kD
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	126961/333932/6536043020/3021
	Gene Symbol	H3C15.H3C14.H3C13 H3-3A.H3-3B H32_HUMAN H33_HUMAN H3C_HUMAN <
	Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Histone H3 (Tri Methyl Lys37)
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.