

## Anti-Acetyl-H2AZ1-Lys14 antibody (STJA0005939)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                |
| Concentration       | 1 mg/mL                                                                                                   |
| Conjugation         | Unconjugated                                                                                              |
| Purification        | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution Range      | WB 1:1000-2000<br>ELISA 1:5000-20000                                                                      |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                   |
| Isotype             | IgG                                                                                                       |
| Molecular Weight    | Observed: 15kD                                                                                            |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.            |

### TARGET INFORMATION

|             |                                                                                           |
|-------------|-------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">3015</a>                                                                      |
| Gene Symbol | <a href="#">H2AZ1</a>                                                                     |
| UniProt ID  | <a href="#">H2AZ_HUMAN</a>                                                                |
| Specificity | This antibody detects endogenous levels of Human, Mouse, Rat Histone H2A.Z (Acetyl Lys14) |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|