

## Anti-Acetyl-H2AC11; H2AC13; H2AC15; H2AC16; H2AC17-Lys96 antibody (STJA0005956)

| GENERAL INFORMATION    |                     |                                                                                                                         |
|------------------------|---------------------|-------------------------------------------------------------------------------------------------------------------------|
|                        | Product Type        | Primary antibodies                                                                                                      |
|                        | Applications        | WB/ELISA                                                                                                                |
|                        | Host / Source       | Rabbit                                                                                                                  |
|                        | Reactivity          | Human/Mouse/Rat                                                                                                         |
| PRODUCT PROPERTIES     |                     |                                                                                                                         |
|                        | Clonality           | Polyclonal                                                                                                              |
|                        | Concentration       | 1 mg/mL                                                                                                                 |
|                        | Conjugation         | Unconjugated                                                                                                            |
|                        | Purification        | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.               |
|                        | Dilution Range      | WB 1:1000-2000<br>ELISA 1:5000-20000                                                                                    |
|                        | Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                                 |
|                        | Isotype             | IgG                                                                                                                     |
|                        | Molecular Weight    | Observed: 15kD                                                                                                          |
|                        | Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                           |
| TARGET INFORMATION     |                     |                                                                                                                         |
|                        | Gene ID             | <a href="#">8329/8330/8332/8336/8969</a>                                                                                |
|                        | Gene Symbol         | <a href="#">H2AC11.H2AC13.H2AC15.H2AC16.H2AC17</a>                                                                      |
|                        | UniProt ID          | <a href="#">H2A1_HUMAN</a>                                                                                              |
|                        | Specificity         | This antibody detects endogenous levels of Human, Mouse, Rat Histone H2A (Acetyl Lys96)                                 |
| ADDITIONAL INFORMATION |                     |                                                                                                                         |
|                        | Note                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |