

## Anti-PPIB antibody [2B10] {Biotin} (STJA0006033)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/WB/IHC          |
| Host / Source | Mouse              |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                 |
| Concentration       | 1 mg/mL                                                                                                                                                    |
| Conjugation         | Biotin                                                                                                                                                     |
| Purification        | The antibody was purified using affinity-chromatography using specific immunogen.                                                                          |
| Dilution Range      | Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: WB 1:1000-2000 IHC 1:50-300 |
| Formulation         | Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.                                                                                       |
| Isotype             | IgG                                                                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                             |

### TARGET INFORMATION

|             |                                                                                                            |
|-------------|------------------------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">5479</a>                                                                                       |
| Gene Symbol | <a href="#">PPIB</a>                                                                                       |
| UniProt ID  | <a href="#">PPIB_HUMAN</a>                                                                                 |
| Specificity | Cyclophilin B Monoclonal Antibody (2B10) Biotin conjugated specially designed for your WB or IHC analysis. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|