

Anti-PPIB antibody [2B10] {Cy3} (STJA0006034)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	IF/WB/IHC
Host / Source	Mouse
Reactivity	Human/Rat/Mouse
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Cy3
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:50-300 IF 1:50-200
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	5479
Gene Symbol	PPIB
UniProt ID	PPIB_HUMAN
Specificity	Cyclophilin B Monoclonal Antibody (2B10) Cy3 Conjugated specially designed for your Immunofluorescence analysis.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.