

## Anti-Flag-Tag antibody [2C5] {AbFluor™ 350} (STJA0006038)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA/IP/IF
Host / Source	Mouse
Reactivity	Species independent

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	AbFluor™ 350
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IF 1:2000
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Specificity	Flag-Tag Monoclonal Antibody (2C5) AbFluor™ 350 Conjugated specially designed for your Immunofluorescence analysis.
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---