

Anti-GAPDH antibody [2B8] {AbFluor™ 405} (STJA0006054)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IF/IHC
Host / Source	Mouse
Reactivity	Human/Rat/Mouse/Monkey/Dog/Chicken/Hamster/Rabbit/Pig/Sheep/Insect/Yeast
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	AbFluor™ 405
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:200 IF 1:200
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2597
Gene Symbol	GAPDH
UniProt ID	G3P_HUMAN
Specificity	GAPDH Monoclonal Antibody (2B8) AbFluor™ 405 Conjugated specially designed for your Immunofluorescence analysis.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.