

Anti-HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B antibody {AbFluor™ 488} (STJA0006093)

GENERAL INFORMATION		
Product Type	Primary antibodies	
Applications	WB	
Host / Source	Mouse	
Reactivity	Zebrafish	
PRODUCT PROPERTIES		
Clonality	Monoclonal	
Concentration	1 mg/mL	
Conjugation	AbFluor™ 488	
Purification	The antibody was purified using affinity-chromatography using specific immunogen.	
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:50-300 IF 1:200-400	
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.	
Isotype	IgG1	
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.	
TARGET INFORMATION		
Gene ID	3020/3021 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968	
Gene Symbol	H3-3A H3-3B H3C1.H3C2.H	
UniProt ID	H33_HUMAN H31_HUMAN H32_HUMAN <	
Specificity	Histone H3 Monoclonal Antibody AbFluor™ 488 Conjugated specially designed for your Immunofluorescence analysis.	
ADDITIONAL INFORMATION		
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.	