

Anti-HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B antibody {Biotin} (STJA0006098)

GENERAL INFORMATION		
Product Type	Primary antibodies	
Applications	WB	
Host / Source	Mouse	
Reactivity	Zebrafish	
PRODUCT PROPERTIES		
Clonality	Monoclonal	
Concentration	1 mg/mL	
Conjugation	Biotin	
Purification	The antibody was purified using affinity-chromatography using specific immunogen.	
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: WB 1:2000-5000 IHC 1:50-300	
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.	
Isotype	IgG1	
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.	
TARGET INFORMATION		
Gene ID	3020/3021 126961/333932/653604 H3-3A H3-3B H3C15 H3C14 H3C13 H33_HUMAN H32_HUMAN H31_HUMAN <	
Specificity	Histone H3 Monoclonal Antibody Biotin conjugated specially designed for your WB or IHC analysis.	
ADDITIONAL INFORMATION		
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.	