

**Anti-  
 HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B  
 antibody {FITC} (STJA0006101)**

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Zebrafish
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are IF (1:250-1:2000) FCM (1:250-1:2000)
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	<a href="#">126961/333932/653604/8350/8351/8352/8353/8354/8355/</a>
Gene Symbol	<a href="#">H3C15.H3C14.H3C13</a>
UniProt ID	<a href="#">H32_HUMAN</a> <a href="#">H31_HUMAN</a> <a href="#">H33_HUMAN</a> <
Specificity	Histone H3 Monoclonal Antibody FITC conjugated specially designed for your WB or IHC analysis.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.