

Anti-LMNB1 antibody [7C11] {HRP} (STJA0006120)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/IP
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: WB 1:2000-5000 IHC 1:50-300
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Molecular Weight	Calculated: 66kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4001
Gene Symbol	LMNB1
UniProt ID	LMNB1_HUMAN
Specificity	Lamin B1 Monoclonal Antibody (7C11) HRP Conjugated, specially designed for your Western blot analysis.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---