

## Anti-MBP-Tag antibody [7G2] {HRP} (STJA0006121)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB                  |
| Host / Source | Mouse               |
| Reactivity    | Species independent |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | HRP  |
| Purification        | The antibody was purified using affinity-chromatography using specific immunogen.  |
| Dilution Range      | Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: WB 1:5000 |
| Formulation         | Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.   |
| Isotype             | IgG  |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.  |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Specificity | MBP-Tag Monoclonal Antibody (7G2) HRP Conjugated, specially designed for your Western blot analysis. |
|-------------|--|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|