

Anti-PCNA antibody [12D10] {AbFluor™ 647} (STJA0006142)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IF/IHC
Host / Source	Mouse
Reactivity	Human/Rat/Mouse
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	AbFluor™ 647
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:200 IF 1:200
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Molecular Weight	Observed: 30-33kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	5111
Gene Symbol	PCNA
UniProt ID	PCNA_HUMAN
Specificity	PCNA Monoclonal Antibody (12D10) AbFluor™ 647 Conjugated specially designed for your Immunofluorescence analysis.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.