

## Anti-V5-Tag antibody [3C8] {FITC} (STJA0006169)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB/IP/IF            |
| Host / Source | Mouse               |
| Reactivity    | Species independent |

### PRODUCT PROPERTIES

|                     |                                                                                                                                          |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                               |
| Concentration       | 1 mg/mL                                                                                                                                  |
| Conjugation         | FITC                                                                                                                                     |
| Purification        | The antibody was purified using affinity-chromatography using specific immunogen.                                                        |
| Dilution Range      | Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IF 1:1000 |
| Formulation         | Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.                                                                     |
| Isotype             | IgG                                                                                                                                      |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                           |

### TARGET INFORMATION

|             |                                                                                                  |
|-------------|--------------------------------------------------------------------------------------------------|
| Specificity | V5-Tag Monoclonal Antibody (3C8) FITC conjugated specially designed for your WB or IHC analysis. |
|-------------|--------------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|