

Anti-VSV-G-Tag antibody (501-511) [8D6] {Cy3} (STJA0006178)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IF
Host / Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Cy3
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IF 1:1000
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	501-511
Specificity	VSV-G-Tag Monoclonal Antibody (8D6) Cy3 Conjugated specially designed for your Immunofluorescence analysis.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---