

Anti-CA9 antibody [PTR2373] (STJA0006249)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was purified by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC-p 1:200-400
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Molecular Weight	Calculated: 49kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	768
Gene Symbol	CA9
UniProt ID	CAH9_HUMAN
Specificity	CA IX monoclonal antibody detects endogenous levels of CA IX protein.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunohistochemistry - Paraffin (IHC-P): Human stomach tissue was stained with Anti-CA IX (PTR2373) Antibody. STJA0006249