

Anti-VDR antibody (150-250) [PTR2542] (STJA0006257)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified using Protein G purification.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS pH7.4 containing 50% Glycerol and 0.03% Proclin 300.
Isotype	IgG1k
Molecular Weight	Calculated: 48kD Observed: 37kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7421
Gene Symbol	VDR
UniProt ID	VDR_HUMAN
Immunogen Region	150-250
Specificity	This antibody detects endogenous levels of Vitamin D Receptor at Human

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Whole cell lysates of Hela were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Vitamin D Receptor (PTR2542) antibody. The HRP-conjugated Goat anti- mouse IgG (H + L) antibody was used to detect the antibody. Lane 1: Hela
 Predicted band size: 48kDa
 Observed band size: 37kDa