

Anti-RRAD antibody (1-100) [PTR2260] (STJA0006261)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol and 0.05% Proclin 300.
Isotype	IgG1k
Molecular Weight	Calculated: 33kD Observed: 33kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6236
Gene Symbol	RRAD
UniProt ID	RAD_HUMAN
Immunogen Region	1-100
Specificity	This antibody detects endogenous levels of RRAD at Human

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Whole cell lysates of A431 were separated by 10% SDS-PAGE, and the membrane was blotted with anti-GSK3β (PTR2553) antibody. The HRP-conjugated Goat anti- mouse IgG (H + L) antibody was used to detect the antibody. Lane 1: A431 Predicted band size: 46kDa Observed band size: 42kDa



Whole cell lysates of A431 were separated by 10% SDS-PAGE, and the membrane was blotted with anti-GSK3β (PTR2553) antibody. The HRP-conjugated Goat anti- mouse IgG (H + L) antibody was used to detect the antibody. Lane 1: A431 Predicted band size: 46kDa Observed band size: 46kDa



Whole cell lysates of HeLa were separated by 10% SDS-PAGE, and the membrane was blotted with anti-GSK3β (PTR2553) antibody. The HRP-conjugated Goat anti- mouse IgG (H + L) antibody was used to detect the antibody. Lane 1: HeLa Predicted band size: 46kDa Observed band size: 46kDa