

## Anti-TLR1 antibody (STJA0006697)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/IF             |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | IHC 1:200-500<br>IF 1:50-200                                                                                          |
| Formulation         | Liquid in PBS pH7.4, containing 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.                                        |
| Isotype             | IgG                                                                                                                   |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION

|             |                                               |
|-------------|-----------------------------------------------|
| Gene ID     | <a href="#">7096</a>                          |
| Gene Symbol | <a href="#">TLR1</a>                          |
| UniProt ID  | <a href="#">TLR1_HUMAN</a>                    |
| Specificity | The antibody detects endogenous TLR1 protein. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Immunohistochemical analysis of paraffin-embedded mouse Spleen Tissue using TLR1 Polyclonal Antibody.